

Ear Spray

Ear cleansing spray that prevents and removes the buildup of earwax



FACT BOX

Ear Spray

Market segments: OTC

Composition: Sterile isotonic seawater solution or alternatively: sterile isotonic saline solution

Product approvals: Registered according to MDD (Medical Device Directive, 93/42/EEC, Class IIa)

Shelf life: 3 years from production date.

Packaging: Standard can sizes are: 50, 100, 125, 150, 200 and 250 ml. Individual carton boxes, with or without leaflets are also available.

Presentations: Both sticker labels and preprinted cans can be used. Labels can be provided with a tamper seal. Preprinted cans are recommended for quantities of 25,000+ per order.

The ear cleansing spray from Aurena Laboratories is a medical device that cleans the ear and removes the buildup of earwax in a safe and convenient way.

Ear Spray help cleanse the ear from wax or debris. The spray also clean ears without harming the ear together with preventing the buildup of more earwax and thus helps to maintain clean and healthy ears. The spray is easy to use and can be sprayed in any angle until the last drop. The formula contains no harsh chemicals and is preservative free.

Since the propellant, compressed air or nitrogen, is not flammable, the spray is classed as a non-flammable aerosol, which offers advantages regarding safety, transportation and storage.

Description:

Remove and prevent the buildup of earwax.
Can be used from 30 months.

Claims:

- Prevention of build up of ear wax
- Removal of ear wax
- Rinsing away ear wax, without pressure to the ear drum
- Helps maintain healthy ears
- Sterile until last drop



AURENA
LABORATORIES

Fjärrviksvägen 22
653 50 Karlstad
Sweden