

Eyewash Spray

Handheld emergency eyewash spray



FACT BOX

Eyewash Spray

Market segments:

OTC, occupational safety

Composition: Sterile 0.9% NaCl solution, Sterile 0,9% NaCl phosphate buffered solution or alternatively sterile Ringer solution.

In addition, actuators for different flow rates and impact pressures are available

Product approvals: Registered according to MDD (Medical Device Directive, 93/42/EEC, Class IIa)

Shelf life: 3 years from production date.

Packaging: Standard can sizes are: 50, 100, 125, 150, 200, 250 and 300 ml. Individual carton boxes, with or without leaflets are also available.

Presentations: Both sticker labels and preprinted cans can be used. Labels can be provided with a tamper seal. Preprinted cans are recommended for quantities of 25,000+ per order.

Accessories: Wall station (for 3 cans) and one can holder

Hand-held saline based eye wash spray for first aid treatment

Aurena Eye Wash is a hand-held sterile, isotonic saline based eye wash spray. The spray is to be used for rinsing and cleansing after exposure to chemically harmful substances or particulate non-penetrating matter or a combination thereof.

The fluid in the device is hypoallergenic, sterilized saline or sterilized Ringer's solution. It may be buffered or nonbuffered.

The device is intended to be used in institutions and workplaces and can be placed in a first aid wall station. It can be used without technical experience or training and has no age limitations. A tamper seal provides highly visible notification if the product has been used or not. Thanks to the separation of product and propellant, the bag-on-valve spray aerosol offers a more constant pressure and flow rate until the last drop, making it possible to spray from any angle.

Since the propellant, compressed air or nitrogen, is not flammable, the spray is classed as a non-flammable aerosol, which offers advantages regarding safety, transportation and storage.

This product contains no preservatives and may be used by all age groups.

Description:

Sterile isotonic saline solution for rinsing/cleansing of the Eye. Flushing of eye can be given constantly for ≥ 15 min with one (250 ml can)

For all ages.

Claims:

- The liquid is suitable for rinsing and cleansing dust and dirt, chemicals, biological hazards and venom from eyes.
- Moisturizing of eye
- Can be used by persons wearing contact lenses
- Hypoallergenic
- Preservative free
- The spray works in any angle.
- Sterile until last drop



Fjärrviksvägen 22
653 50 Karlstad
Sweden